Better products • Better care ♥

# ForaCare<sup>®</sup> GD40

Blood Glucose Monitoring System

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# Dear ForaCare GD40 System Owner:

Thank you for purchasing the **ForaCare GD40** Blood Glucose Monitoring System. This manual provides important information to help you to use the system properly. Before using this product, please read and follow the instructions in this Owner's Manual and any package inserts for test strips and control solution.

Regular monitoring of your blood glucose levels can help you and your doctor gain better control of your diabetes. Due to its compact size and easy operation, you can use the **ForaCare GD40** Blood Glucose Monitoring System to easily monitor your blood glucose levels.

If you have questions, just call ForaCare Customer Care Service Center toll-free at 1-866-469-2632. We are here for you 24 hours a day, 7 days a week. You can also visit [www.foracare.com](http://www.foracare.com).

# **IMPORTANT SAFETY PRECAUTIONS READ BEFORE USE**

- The meter and lancing device are for single patient use. Do not share them with anyone including other family members! Do not use on multiple patients!
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.

For more information, please visit

1. "FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication" (2010)

<http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>

2. "CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens" (2010)

<http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>

1. Use this device **ONLY** for the intended use described in this manual.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or if it is damaged.
4. Do **NOT** use the equipment in places where aerosol sprays are being used or where oxygen is being administered.
5. Do **NOT** under any circumstances use the device on neonates.
6. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only.
7. Before using this device to test blood glucose, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
8. Keep the device and testing equipment away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
9. Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
10. Do **NOT** use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.

**KEEP THESE INSTRUCTIONS IN A SAFE PLACE**

# TABLE OF CONTENTS

<b>BEFORE YOU BEGIN</b>	<b>08</b>
Important Information	08
Intended Use	09
Test Principle	09
Meter Overview	10
Display Screen	11
Test Strip	12
Setting the Meter	13
<b>THE FOUR MEASURING MODES</b>	<b>17</b>
<b>BEFORE TESTING</b>	<b>18</b>
Calibration	18
How to Code Your Meter	18
Checking the Code Number	19
Control Solution Testing	20
Performing a Control Solution Test	20
<b>TESTING WITH BLOOD SAMPLE</b>	<b>23</b>
Preparing the Lancing Device for Blood Testing	23
Preparing the Puncture Site	24
Performing a Blood Glucose Test	25
<b>METER MEMORY</b>	<b>28</b>
<b>DOWNLOADING RESULTS ONTO A COMPUTER</b>	<b>30</b>

<b>MAINTENANCE</b>	<b>32</b>
Battery	32
Caring for Your Meter	34
Meter Storage	37
Caring for Your Test Strips	37
Important Control Solution Information	38
<b>SYSTEM TROUBLESHOOTING</b>	<b>39</b>
Result Readings	39
Error Messages	40
Troubleshooting	41
<b>DETAILED INFORMATION</b>	<b>42</b>
Reference Values	42
<b>SPECIFICATIONS</b>	<b>43</b>
<b>FEDERAL COMMUNICATIONS COMMISSION (FCC) STATEMENT</b>	<b>44</b>

## BEFORE YOU BEGIN

### Important Information

---

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than usual, and you do not have any symptoms of illness, first repeat the test. If you have symptoms or continue to get results which are higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole capillary blood samples from the finger to test your blood glucose. Use only heparin (sodium heparin) for anticoagulation of whole blood. Using other substances will lead to incorrect results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all the instructions given in this owner's manual, contact your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis. Please consult the healthcare professional before use.

## Intended Use

---

This system is intended for use outside the body (*in vitro* diagnostic use) by people with diabetes at home and by health care professionals in clinical settings as an aid in monitoring the effectiveness of diabetes control. It is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips. It should not be used for the diagnosis of or screening for diabetes mellitus.

It should not be used for the diagnosis of or screening for diabetes, or testing on neonates.

## Test Principle


---

Your system measures the amount of sugar (glucose) in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

This system uses enzyme of FAD-dependent glucose dehydrogenase (FAD-GDH) from *E. coli*.

## Meter Overview

---

**1 Test Slot**

Insert test strip here to turn the meter on for testing.

**2 M Button**

Enter the meter memory and silence a reminder alarm.

**3 Display Screen**

**4 Test Strip Ejector**

Eject the used strip by pushing up this button.

**5 Data Port** (for GD40a only)

Download test results with a cable connection.


**6 SET Button**

Enter SET mode and confirm the meter settings.

**7 Battery Compartment**

## Display Screen

---

1 Blood drop symbol

2 Code

3 Test Result

4 Ketone warning

5 Error Warning

6 Memory symbol

7 Measurement modes

8 Date

9 Test strip symbol

10 Low battery symbol

11 Measurement unit


12 Reminder alarms

13 Day average

14 Time

## Test Strip

---

### ATTENTION:

**The front side of test strip should face up when inserting test strip.**

Test results might be wrong if the contact bar is not fully inserted into the test slot.

### NOTE:

The **ForaCare GD40** meter should only be used with **ForaCare GD40** Test Strips. Using other test strips with this meter can produce inaccurate results.


## Setting the Meter

---

Before using your meter for the first time or if you change the meter battery, you should check and update these settings. Make sure you complete the steps below and have your desired settings saved.

### ► Entering the Setting Mode


Start with the meter off (no test strip inserted). Press SET to turn on the meter.

#### Step 1

#### Setting the date

With the year flashing, press **M** until the correct year appears. Press SET button.

With the month flashing, press **M** until the correct month appears. Press SET button.


With the day flashing, press **M** until the correct day appears. Press SET button.

## Step 2

### Setting the time format

Press and release **M** to select the desired time format — 12h or 24h.  
Press SET button.

## Step 3

### Setting the time

With the hour flashing, press **M** until the correct hour appears. Press SET button.

With the minute flashing, press **M** until the correct minute appears. Press SET button.


## Step 4

### Deleting the memory

With "dEL" and a flashing "**M**" symbol on the display, press **M** and select "No" to keep the results in memory then press SET to skip.


To delete all the results, press **M** and select “Yes” to delete the results in memory, then press SET to confirm. “OK” and “**M**” are displayed on the meter, which indicates that all data stored is deleted.


### Step 5

#### Setting the reminder alarm

You may set up any or all of the reminder alarms (1-4). The meter displays “On” or “OFF” and <sup>1</sup>🔔, press **M** to turn on or turn off to set the first reminder alarm.


Press **M** to select “On”, then press SET to set the hour. When the hour is flashing, press **M** to add an hour. Press SET to confirm and go to minutes, press **M** to add one minute. Hold **M** longer to add faster. Press SET to confirm and go to the next alarm setting.

If you do not want to set an alarm, press SET to skip this step.  
If you want to turn off an alarm, find the alarm number by pressing SET in the setting mode, press **M** to change from “ON” to “OFF”.

At the time of your alarm, the meter will beep and automatically turn on. You can press **M** to silence the alarm and insert a test strip to begin testing.

If you do not press **M**, the meter will beep for 2 minutes then switch off.

If you do not want to test at this time, press **M** to switch off the meter.

**Congratulations! You have completed all settings!**

#### NOTE

- These parameters can **ONLY** be changed in the setting mode.
- If the meter is idle for 3 minutes during the setting mode, it will switch off automatically.

## THE FOUR MEASURING MODES

The meter provides you with four modes for measuring, General, AC, PC and QC.

Modes	Use when
General (displays as "Gen")	any time of day without regard to time since last meal
<b>AC</b>	no food intake for at least 8 hours
<b>PC</b>	2 hours after a meal
<b>QC</b>	testing with the control solution

You can switch between each mode by:


### Step 1

Start with the meter switched off. Insert a test strip to turn on the meter, the screen will display flashing blood drop and "Gen".

### Step 2

Press **M** to switch between General, AC, PC and QC mode.

## BEFORE TESTING

### Calibration

---


You must calibrate the meter every time you begin to use a new vial of test strips by setting the meter with the correct code. Test results may be inaccurate if the code number displayed on the monitor does not match the number printed on the strip vial.

### How to Code Your Meter

---

#### Step 1

Insert the code strip when the meter is off. Wait until the code number appears on the display.


#### NOTE

Make sure the code number on display, code strip, and test strip vial are the same.

## Step 2


Remove the code strip, the display will show “OK”.  
This tells you that the meter has finished coding and is ready for blood glucose testing.

## Checking the Code Number

---

You need to make sure that the code number displayed on the meter matches the number on the test strip vial before you proceed. If it matches, you can proceed with your test. If the codes do not match, please stop testing and repeat the calibration procedure. If the problem persists, contact Customer Service for help.

### Warning

It is important to make sure that the LCD displayed code is the same as the code on the test strip vial before testing. Failure to do so will get inaccurate results.

## Control Solution Testing

---

Our Control Solution contains a known amount of glucose that reacts with test strips and is used to ensure your meter and test strips are working together correctly.

### Do a control solution test when:

- you first receive the monitor
- at least once a week to routinely check the monitor and test strips
- you begin using a new vial of test strips
- you suspect the monitor or test strips are not working properly
- your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate
- practicing the testing process, or
- you have dropped or think you may have damaged the monitor.


## Performing a Control Solution Test

---

### Step 1

#### Insert the test strip to turn on the meter


Insert the test strip into the meter. There will be a strip light indicating the test strip is inserted into the test slot. Wait for the meter to display the test strip and blood drop symbols.

## Step 2

### Press **M** to Mark this Test as a Control Solution Test.

With “QC” displayed, the meter will store your test result in memory under “QC”. If you press the **M** again, the “QC” will disappear and this test is no longer a control solution test.

## Warning


When doing the control solution test, you have to mark it so that the test result will not mix with the blood glucose stored in the memory.

## Step 3

### Apply Control Solution.

Shake the control solution vial thoroughly before use. Squeeze out a drop and wipe it off, then squeeze another drop and place it on the tip of the vial cap.

Hold the meter to move the absorbent hole of test strip to touch the drop. Once the confirmation window fills completely, the meter will begin counting down.


To avoid contaminating the control solution, do not directly apply control solution onto a strip.

#### Step 4

#### Read and Compare the Result

After counting down to 0, the control solution test result will appear on the display. Compare this result with the range printed on the test strip vial and it should fall within this range. If not, please read the instructions again and repeat the control solution test.

#### Out-of-range results

If you continue to have test results fall outside the range printed on the test strip vial, the meter and strips may not be working properly. Do **NOT** test your blood. Contact the local customer service or place of purchase for help.

#### NOTE

- The control solution range printed on the test strip vial is for control solution use only. It is not a recommended range for your blood glucose level.
- See the **MAINTENANCE** section for important information about your control solutions.

## TESTING WITH BLOOD SAMPLE

### ► Preparing the Lancing Device for Blood Testing

Please follow the instructions in the lancing device insert for collecting a blood sample.

#### WARNING

To reduce the chance of infection:

- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.
- Wash and dry your hands thoroughly after handling the meter, lancing device and test strips to prevent infection. For more information, please refer to the “Cleaning and Disinfection” section.
- If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be decontaminated prior to use by the second person.

**Sharing the lancing device and lancets may increase the risk of contracting infectious diseases. Lancing device must not be used for more than one person.**

## Preparing the Puncture Site

---


Stimulating blood perfusion by rubbing the puncture site before blood extraction has a significant influence on the glucose value obtained. Blood from a site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site was rubbed prior to blood extraction, the difference was significantly reduced.

### **Please follow the suggestions below before obtaining a drop of blood:**

- Wash and dry your hands before starting.
- Select the puncture site at fingertips.
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry**.
- Rub the puncture site for about 20 seconds before penetration.

### **► Fingertip testing**

Press the lancing device's tip firmly against the lower side of your fingertip. Press the release to prick your finger, then a click indicates that the puncture is complete.

### **NOTE**

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- It is recommended to discard the first drop of blood as it might contain tissue fluid, which may affect the test result.


## Performing a blood glucose Test

---

### Step 1

#### Insert the test strip to turn on the meter

There will be a strip light indicating the test strip is inserted into the test slot. Wait for the meter to display the "  " and "  " symbols.

### Step 2


#### Select the appropriate measuring mode by pressing **M**.


For selecting the measurement mode, please refer to the "FOUR MEASURING MODES".

### Step 3

#### Obtaining a blood sample

Use the pre-set lancing device to puncture the desired site. After penetration, discard the first drop of blood with a clean cotton swab. Gently squeeze the punctured area to obtain another drop of blood. Be careful **NOT** to smear the blood sample.



The volume of blood sample must be at least **1.1** microliter ( $\mu\text{L}$ ) of volume. (  actual size).

#### Step 4

#### Apply the sample

Hold the meter with the strip pointing downward and the absorbent hole of the test strip touching the blood drop. Blood will be drawn in to the strip. Once the confirmation window is completely filled, meter begins counting down. Do not remove the strip from the blood drop until the meter beeps.

#### NOTE

- Do not press the site against the test strip or try to smear the blood.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to start a new test.
- The confirmation window should be filled with blood before the meter begins to count down. **NEVER** try to add more blood to the test strip after your drop of blood has moved away. **Discard the used test strip and retest with a new one.**
- If you have trouble filling the confirmation window, please contact your health care professional or the local customer service for assistance.

### Step 5

#### Read the result


The result of your blood glucose test will appear after the meter counts down to 0. This blood glucose result will automatically be stored in the memory.

### Step 6

#### Eject the used test strip

To eject the test strip, point the strip into a disposal container for contaminated items. Slide the Test Strip Ejector forward. The meter will automatically turn itself off after the test strip is ejected.

**Always follow the instructions in the lancing device insert when removing the lancet.**

### WARNING

The used lancet and test strip may be biohazards. Please discard them carefully according to your local regulations.

Wash your hands thoroughly with soap and water after handling the meter, lancing device and test strips to avoid contamination. For more information, please refer to the “Cleaning and Disinfection” section.

## METER MEMORY


The meter stores the 1000 most recent blood glucose test results along with respective dates and times in its memory. To enter the meter memory, **start with the meter switched off.**

### ► Reviewing Test Results

#### Step 1


Press and release **M**.

“**M**” will appear on the display. Press **M** again, and the first reading you see is the last blood glucose result along with date, time and the measuring mode.

#### Step 2

Press **M** to recall the test results stored in the meter each time you press.

#### Step 3

Exit the meter memory.


After the last test result, press **M** again and the meter will be turned off.

## ► Reviewing Blood Glucose Day Average Results

### Step 1

#### Press and release **M**

When **M** appears on the display, keep pressing **M** for 3 seconds until the flashing “**DAY AVG**” appears. Release **M** and then your 7-day average result measured in general mode will appear on the display.

### Step 2

Press **M** to review 14-, 21-, 28-, 60- and 90- day average results stored in each measuring mode in the order of Gen, AC, then PC.


### Step 3

#### Exit the meter memory

Keep pressing **M** and the meter will turn off after displaying the last test result.

### NOTE

- Any time you wish to exit the memory, keep pressing **M** for 5 seconds or leave it without any action for 3 minutes. The meter will switch off automatically.
- Control solution results are **NOT** included in the day average.
- If using the meter for the first time, “—” displays when you recall the test results or review the average result. It indicates that there is no test result in the memory.

## DOWNLOADING RESULTS ONTO A COMPUTER

### ► Data transmission via cable (For Foracare GD40a)

You can use the meter with the USB cable and the Health Care Software System to view test results on your personal computer. To learn more about the Health Care Software System or to obtain the USB cable separately, please contact local customer services or the place of purchase for assistance.

#### Step 1


##### Obtaining the required cable and installing the software

To download the Health Care Software System, please visit the Foracare Inc. website: <http://www.foracare.com>.

#### Step 2

##### Connecting to a personal computer

Connect the cable to the USB port on your computer. With the meter switched off, connect the other end of the USB cable to the meter data port. "USb" will appear on the meter display, indicating that the meter is in communication mode.

#### Step 3

##### Data transmission

To transmit data, follow the instructions provided with the software. Results will be transmitted with date and time. Remove the cable and the meter will automatically switch off.

### WARNING

While the meter is connecting to the PC, it will be unable to perform a blood glucose test.

### ► **Data transmission via Bluetooth (For Foracare GD40b)**

You can use your monitor and Health Care Software System to view your test results on your personal computer. To learn more about the Health Care Software System, please contact local customer services or the place of purchase for assistance.

#### **Step 1**


##### **Installing the software**

To download the Health Care Software System, please visit the Foracare Inc. website: <http://www.foracare.com>.

#### **Step 2**

##### **Connecting to a personal computer**

The meter will automatically initiate searching for the compatible Bluetooth receiver and start the data transmission after the measurement, the calibration, exit the meter memory or finishing the setting. The BT indicator will light up in blue and start blinking.

#### **Step 3**

##### **Data transmission**

When the BT indicator is blinking, this indicates that the meter is searching for pairing. When the BT indicator shows a steady blue light, the meter has connected to the Bluetooth receiver and the transmission is going to begin.

The BT indicator will blink for 3 minutes. If the compatible Bluetooth receiver was not found, then the indicator would be off automatically afterwards.

#### **NOTE**

- Make sure your PC has Bluetooth before transmitting the data and the meter is within the receiving range.
- While the meter is connected to the PC, it is unable to test for blood glucose.

## MAINTENANCE

### Battery

---

Your meter comes with two 1.5V AAA size alkaline batteries.


#### ► Low Battery Signal

The meter will display one of the messages below to alert you when the meter power is getting low.

##### Step 1

The “” symbol appears along with display messages:


The meter is functional and the result remains accurate, but it is time to change the batteries.

##### Step 2

The  symbol appears with **E-b** and  **LOW** :

The power is not enough to do a test. Please change the batteries immediately.

## ► Replacing the Battery

To replace the battery, make sure that the meter is turned off.

### Step 1


Press the edge of the battery cover and lift it up to remove.

### Step 2

Remove the old batteries and replace with two 1.5V AAA alkaline batteries.

### Step 3

Close the battery cover. If the batteries are inserted correctly, you will hear a "beep" afterwards.

### NOTE

- Replacing the batteries does not affect the test results stored in the memory.
- As with all small batteries, these batteries should be kept away from small children. If swallowed, promptly seek medical assistance.
- Batteries might leak chemicals if unused for a long time. Remove the batteries if you are not going to use the device for an extended period (i.e., 3 months or more).
- Properly dispose of the batteries according to your local environmental regulations.
- After replacing the batteries, the meter will enter the setting mode.

## Caring for Your Meter

---

To avoid the meter and test strips attracting dirt, dust or other contaminants, please wash and dry your hands thoroughly before use.

### ► **Why the cleaning and disinfection should be performed**

Cleaning and disinfection are different. Cleaning is the process of removing dirt (e.g. food debris, grease, dust), disinfection is the process of killing germs (e.g. bacteria and viruses).

### ► **When to clean and disinfect the meter**

Clean the meter when you see any dirt on it. You should disinfect the meter at least once a week to prevent infection.

### ► **How to clean and disinfect the meter**

The meter must be cleaned prior to the disinfection. Use one disinfecting wipe to clean exposed surfaces of the meter thoroughly and remove any visible dirt or blood or any other body fluid with the wipe. Use a second wipe to disinfect the meter. Do **NOT** use organic solvents to clean the meter.

We recommend for meter cleaning and disinfection you should use the disinfecting wipes/towelettes from below. The following product with isopropyl alcohol as the active ingredient has been shown to be safe for use with the ForaCare GD40 Blood Glucose Monitoring System.

## ► Disinfecting Procedures

1. Take out one disinfecting wipe from the package and squeeze out any excess liquid in order to prevent damage to the meter.

2. Wipe all meter's exterior surface display and buttons. Hold the meter with the test strip slot pointing down and wipe the area around the test slot but be careful not to allow excess liquid to get inside. Keep meter wet with disinfection solution contained in the wipe for a minimum of 2 minutes for Micro-Kill+™ wipes. Follow the instructions on the package label of disinfecting wipe.

3. Remove the wipe. Allow the meter surface to dry completely.

4. Discard the used wipes and never reuse them.

Improper system cleaning and disinfection may result in meter malfunction. If you have a question, please contact local customer service at 1-866-469-2632 for assistance.

This device has been validated to withstand up to 10,000 cleaning and disinfection cycles using the recommended disinfecting wipe/towelette. The tested number of cycles is estimated by 5 cleaning and disinfection cycles per day over 5 years, the expected life of the device. The meter should be replaced after the validated number of

cleaning and disinfection cycles or the warranty period, whichever comes first.

Stop using the meter if you see any signs of deterioration. For example:

- meter cannot be turned on,
- LCD display cracks or becomes cloudy,
- buttons no longer function,
- meter outer casing cracks,
- data can not be transmitted to PC,
- color or paint/printing on housing is abnormal, or
- scratches or abrasions on meter are higher than acceptable.

Please contact the customer service for a replacement meter if any of the signs of deterioration are noticed.

#### **NOTE**

- Do **NOT** clean and disinfect the meter while performing tests.
- If the meter is being operated by a second person, the meter and lancing device should be decontaminated prior to use by the second person.
- Do **NOT** allow cleaning and disinfecting solution to get in the test slot, battery compartment, or strip-ejection button.
- If you do get moisture in the test strip slot, wipe it away with a corner of tissue.
- Always dry the meter thoroughly before using it
- Do not spray the meter directly with cleaning solutions especially those containing water (i.e. soapy water), as this could cause the solution to enter the case inside and damage the electronic components or circuitry

## Meter Storage

---

- Storage conditions: -4 °F to 140 °F (-20 °C to 60 °C), below 95% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

## Caring for Your Test Strips

---

- Storage conditions: 35.6 °F to 89.6 °F (2 °C to 32 °C) below 85% relative humidity. Do not freeze.
- Store your test strips in their original vial only. Do not transfer to another container.
- Store test strip packages in a cool dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the vial label when you first opened it. Discard remaining test strips after 3 months.
- Do not use test strips beyond the expiry date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

For further information, please refer to the test strip package insert.

## Important Control Solution Information



---

- Use only **ForaCare** control solutions with your meter.
- Do not use the control solution beyond the expiry date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.
- It is recommended that the control solution test be done at room temperature 68° F to 77° F (20° C to 25° C). Make sure your control solution, meter, and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 35.6° F to 86° F (2° C to 30° C). Do **NOT** freeze.


## SYSTEM TROUBLESHOOTING

If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local customer service. Do not attempt to repair by yourself and never try to disassemble the meter under any circumstances.

### Result Readings

MESSAGE	WHEN GLUCOSE
	Appears when your result is below measurement limit, which is less than 20 mg/dL (1.1 mmol/L).
	Appears when your result is equal to or higher than 240 mg/dL (13.3 mmol/L). This indicates the possibility of ketone accumulation for type 1 diabetes. Please seek medical assistance immediately.
	Appears when your result is higher than the limit of measurement, which is higher than 600 mg/dL (33.3mmol/L).

## Error Messages

MESSAGE	WHAT IT MEAN	ACTION
 <p><b>E-b</b> ⚠️ <small>LOW</small></p>	Appears when the batteries cannot provide enough power for a test.	Replace the batteries immediately.
 <p><b>E-U</b> ⚠️</p>	Appears when a used test strip is inserted.	Repeat the test with a new test strip.
 <p><b>E-E</b> ⚠️</p>	Problem in operation.	Repeat the test with a new test strip. If the meter still does not work, please contact the customer service for assistance.
 <p><b>E-F</b> ⚠️</p>	Appears when test strip is removed while counting down, or insufficient blood volume.	Repeat the test with a new test strip.
 <p><small>Code</small> <b>E-C</b> ⚠️</p>	Appears when the wrong code strip is inserted or other coding errors.	Make sure to insert the right code strip for the meter. Check if the code numbers on the code strip and the strip vial label are the same.
 <p><b>E-2</b> ⚠️</p>	Expired code strip.	Repeat the test with a new code strip.
 <p><b>E-t</b> ⚠️ <small>LOW</small></p>	Appears when ambient temperature is below system operation range.	System operation range is 50°F to 104°F (10°C to 40°C). Repeat the test after the meter and test strip are in the above temperature range.
 <p><b>E-t</b> ⚠️ <small>HIGH</small></p>	Appears when ambient temperature is above system operation range.	

## Troubleshooting

---

### 1. If the meter does not display a message after inserting a test strip:

POSSIBLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
Defective meter or test strips.	Please contact customer services.

### 2. If the test does not start after applying the sample:

POSSIBLE CAUSE	WHAT TO DO
Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic switch-off (2 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when flashing “  ” appears on the display.
Defective meter.	Please contact customer services.

### 3. If the control solution testing result is out of range.

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read instructions thoroughly and repeat the test again.
Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiry date of the control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should be at room temperature 68°F to 77°F (20°C to 25°C) before testing.
Defective test strip.	Repeat the test with a new test strip.
Meter malfunction.	Please contact customer services.

## DETAILED INFORMATION

### Reference Values

---

The meter provides you with plasma equivalent results.

Time of day	Normal plasma glucose range for people <b>without</b> diabetes (mg/dL)
Fasting and before meal	<100 mg/dL (5.6 mmol/L)
2 hours after meals	<140 mg/dL (7.8 mmol/L)

**Source:** American Diabetes Association (2010). Clinical Practice Recommendations. Diabetes Care, 33 (Supplement 1): S1-S100.

**Please consult your doctor to determine a target range that works best for you.**

## SPECIFICATIONS

**Model No.:** TD-4272

**Dimension & Weight:** 110 (L) x 57 (W) x 25 (H) mm, 71 g

**Power source:** two 1.5V AAA alkaline batteries

**Display:** LCD

**Memory:** 1000 measurement results with respective date and time

**External output:** USB Cable (ForaCare GD40a) or Bluetooth (ForaCare GD40b)

Auto electrode inserting detection

Auto sample loading detection

Auto reaction time count-down

Auto shutdown after 3 minutes of idleness

Temperature Warning

**Operating condition:** 50°F to 104°F (10°C to 40°C), 10% to 85% R.H.  
(no condensing)

**Storage/Transportation condition:** -4°F to 140°F (-20°C to 60°C), 95% R.H.

**Measurement units:** mg/dL

**Measurement range:** 20 to 600mg/dL (1.1 to 33.3mmol/L)

This device has been tested to meet the electrical and safety requirements of:  
IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1, IEC/EN 61326-2-6.

## FEDERAL COMMUNICATIONS COMMISSION (FCC) STATEMENT

15.21

**You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.**

15.105(b)

### **Federal Communications Commission (FCC) Statement**

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

**This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:**

1. This device may not cause harmful interference and
2. This device must accept any interference received, including interference that may cause undesired operation of the device.

### **FCC RF Radiation Exposure Statement:**

1. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.
2. This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 centimeters between the radiator radiation source and your body.

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Fora Care Inc.